

## New Lot of PT Reagent Worksheet

Date(s) of Testing: \_\_\_\_\_ Date of Initial Use: \_\_\_\_\_

Reagent: \_\_\_\_\_ Lot #: \_\_\_\_\_ Expiration Date: \_\_\_\_\_

New ISI Value: \_\_\_\_\_ Date ISI Programmed in Analyzer: \_\_\_\_\_ Tech: \_\_\_\_\_

### Establishment of Normal Patient Mean w/ New Reagent Lot

Normal Patient PT Results					New Patient Mean:  Analyzer Programmed: _____ Date/Tech: _____
1	5	9	13	17	
2	6	10	14	18	
3	7	11	15	19	
4	8	12	16	20	

### Comparison of Patients' PT/INR Results

Patients in Therapeutic PT/INR Range				Patients w/ Abnormal PT/INR Results			
Old Lot Reagent		New Lot Reagent		Old Lot Reagent		New Lot Reagent	
PT	INR	PT	INR	PT	INR	PT	INR
1				1			
2				2			
3				3			
4				4			
5				5			

### Crossover Quality Control (QC) PT Study w/ Old & New Reagent Lots

Control Level 1 Range:				Control Level 2 Range:			
Old Lot	New Lot	Old Lot	New Lot	Old Lot	New Lot	Old Lot	New Lot
1		6		1		6	
2		7		2		7	
3		8		3		8	
4		9		4		9	
5		10		5		10	

#### Manual INR Check

INR = (Patient PT / Mean Normal Range PT)<sup>ISI</sup>

Patient ID #	Analyzer Results	Manual Calculation
1		
2		
3		

#### Platelet Poor Plasma Check

Centrifuge RPM/Time: \_\_\_\_\_

Platelet Count (less than 10,000)
1
2
3
4
5

Supervisor's Review: \_\_\_\_\_ Date: \_\_\_\_\_